

Reveal Your Beauty

Opeau Rose

User manual




Table of contents

Preface	01-02
Product Notice	01-02
Included In Your Package	03-04
Usage Precautions	05-06
Safety and Operations	07-08
Using the Opeau Rose	09-10
Skin Sensitivity Testing	11-12
Choose Your Energy Level	11-12
Treatment Overview	12-13
Warranty	13-14

Warranty Card

Product Name: IPL beauty hair-removal device for home use	Model: Opeau Rose T4
Purchase Date:	
Standard: GB4706.1-2005	

This product has been tested and approved and complies with the standard safety and quality requirements.

►Preface

Thank you for choosing the Opeau Rose IPL beauty device. Opeau is the leader in IPL with adjustable power levels, skin sensing technology to prevent accidental flashing, gliding technology to quickly cover larger areas and an industry leading "ice cooling" feature that calms and protects your skin from any heat build up.


Please read through the user manual carefully before using the device.

We hope you enjoy using your Opeau Rose. If you have any problems using the device, please email **support@opeau.com**.

1. Product Notice

- This user manual will guide you through using your Opeau IPL beauty device.
- Please use the device according to correct operational steps outlined in this manual.


2. Included In Your Package




1 IPL handset



1 pair of goggles



1 Adapter




User manual

Model No.: Opeau T4

- Input Power: 48W
- Voltage Rating: 100V-240V
- Rated Frequency Range: 50Hz/60Hz
- Storage Temperature: 0°C-45°C (32°F-113°F)
- Working Temperature: 5°C-35°C (41°F-95°F)
- Working Humidity: 25%-75%

(Figure 1)



3. Usage Precautions

-  Do not use if you are under the age of 18.
-  Do not use this device on tanned skin as it may result in temporary or even permanent whitening or darkening of the skin.
-  Do not use if you are pregnant or breast-feeding. Hormonal changes occur during pregnancy and breast-feeding and could increase the risk of skin injury or skin sensitivity.
-  Do not use if you have epilepsy or light-induced seizures or migraines. The energy flashes may cause severe headaches and migraines.
-  Do not use on tattoos, dark brown or black spots (such as dark freckles, birthmarks, moles or warts). Tattoos or dark spots may absorb too much light energy which could cause skin swelling, temporary or permanent blisters, discoloration of the skin or whitening of the skin.


-  Do not use the device on sensitive areas of your body including nipples, genitalia, or around the anus. Use in these areas could cause skin swelling, blisters, or temporary skin discoloration or whitening of the skin.
-  Do not use if you have any active skin condition in the treatment area, including acne, sores, psoriasis, eczema, burns, herpes simplex, wounds or a rash. These conditions can make the skin more sensitive to the light treatment and could make the treatment uncomfortable or cause skin swelling, blisters or infection. Please consult with your doctor before using.
-  Do not use the device if you have a history of keloid scar formation, a known sensitivity to light (photosensitivity) or are taking medication which makes the skin more sensitive to light, including non-steroidal anti-inflammatory agents (e.g. aspirin, ibuprofen, acetaminophen), tetracycline, phenothiazine, thiazide, diuretics, sulfonylurea, sulfonamides, DTIC, fluorouracil, vinblastine, griseofulvin, Alpha-Hydroxy Acids (AHAs), Beta-Hydroxy Acids (BHAs), Retin-A, Accutane and/or topical retinoids.
-  Do not use if you have a disease related to photosensitivity, such as porphyria, polymorphic light eruption, solar urticaria or lupus. These conditions make the skin more sensitive to the light treatment and could make the treatment uncomfortable or cause skin swelling, blisters or infection.

4. Safety and Operations


- Wear the goggles provided while operating the device to protect your eyes.
- Your skin must be clean and free from cosmetics and cream.
- As outlined in the sensitivity testing section, always start with a lower intensity level and raise this gradually.
- Use a soft, dry, alcohol-free cloth to clean the device.
- Store the device in a cool, dry and well-ventilated environment.
- Keep it away from children.
- Do not stick any metallic objects into the IPL heating element to avoid risk of electric shock.
- Do not try to disassemble and repair the device yourself.
- If you experience a problem with this device during the warranty period, please contact **support@opeau.com** for assistance.

5. Using the Opeau Rose

- Clean the cartridge light output window with a cotton swab. This is done in order to remove hairs and dust from the treatment window.
- Connect the power into the base of the handle on your Opeau IPL device. The status indicator on the top of the device will flash slowly.
- To turn the unit on, hold down the power button (refer to figure 1) for 2 seconds. The internal fan will start. When the LED light on top of the unit stops flashing, it indicates that the device is ready. The display will show the remaining flash count, the intensity level and whether the ice cooling feature is active (refer to figure 2).
- The Opeau Rose comes with an advanced skin sensing technology. The green status indicator will only light up when the device head is in contact with the skin. This unit can only be operated when the green indicator is visible.



- By default the Opeau Rose T4 will start at the lowest energy level setting 1. Press the power button to adjust the energy level between 1 and 5.
- Hold the device firmly against the treatment area (refer to figure 3) and press the operation trigger. Move along the skin after each flash. For larger body areas use gliding mode. Hold the trigger button for 3 seconds, the energy level will flash on and off. Slowly move the device along the area to be treated, the device will continue to flash. To cancel gliding mode, press the trigger button.
- The device will power off automatically if not used for 5 minutes.



6.Skin Sensitivity Testing

If you are using the system for the first time or on a new body area, choose a small area close to the area you intend to treat, for the purpose of performing a skin sensitivity test. Use level one (low) energy level. Follow the treatment instructions below but make sure to emit no more than 2-3 pulses. After the skin test, wait 60 minutes and check your skin for any reaction. If after 60 minutes the skin appears normal and there is no redness or other skin reaction, then you are ready to perform a full treatment.

7.Choose Your Energy Level

Choose an energy according to your skin tolerance and increase the level from 1 to 5 gradually. Do not flash multiple times in the same area. Doing so may cause too much heat and may result in skin swelling, blisters, or temporary or permanent changes in the color of your skin. After each flash, move the treatment window slowly along.

8.Treatment Overview

- This treatment is not suitable for all skin types. Please ensure that you meet the criteria outlined in figure 4 below. It is most effective on light to medium skin tones with hair ranging from natural blonde to dark brown or black. The device is not effective on very blonde, red, grey or white hair where the smaller amount of melanin would not absorb the light energy.

(Figure 4)

Skin and natural hair color comparison table

Hair color \ Skin color	White	Ivory tint	Beige	Light brown	Brown	Dark brown and black
Black	✓	✓	✓	✓	✓	✗
Dark brown	✓	✓	✓	✓	✓	✗
Brown	✓	✓	✓	✓	✓	✗
Golden	✓	✓	✓	✓	✓	✗
Red	✗	✗	✗	✗	✗	✗
White	✗	✗	✗	✗	✗	✗
Gray	✗	✗	✗	✗	✗	✗

✓ Applicable ✗ Not applicable

- For optimal results, shave and clean the targeted area before beginning the treatment.
- Wear the goggles provided to avoid IPL light leakage that may affect your eyes.
- To choose a suitable energy level, start with a low level and gradually work upwards.
- Do not repeatedly treat the same area.
- Apply moisturizing cream after treatment.
- Use a cotton swab to clean the treatment window after usage.
- When finished, place the device and accessories into its packaging and store it in a safe dry place.
- **Treatment Regime:** Do one treatment a week for 3 weeks. For the subsequent 5 treatments, do one treatment every 2-3 weeks. You can then adjust to your needs (i.e. every 1-2 months). The result is permanent visible hair removal.

9. Warranty

1) Warranty period

- The warranty is 1 year from the date of purchase and must have a valid remaining lamp count of more than 30%. This warranty does not apply to commercial usage.

2) Warranty scope

- You may not make a claim under warranty unless the defect is due to faulty, defective parts or workmanship. If the unit is damaged by accident, misuse or abuse, usage under abnormal conditions, fire, earthquake, flood, disassembly, abnormal wear and tear, power surges, damage or incorrect power supply, improper storage, incorrect or inappropriate operations and failure to comply with the instructions supplied with the appliance, then the warranty is void.
- Within the warranty period, Opeau will replace the handset with a new unit subject to the conditions outlined above.
- When seeking a warranty claim, please contact **support@opeau.com** with a valid proof of purchase.